



## STACKABLE 200µL TIP BOX BPIC-707

## STACKABLE 200?L TIP BOX BPIC-707

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## BPIC-707 STACKABLE 200?L TIP BOX



Made from medical grade Polypropylene  
Slippy surface with excellent design  
DNase-free  
Rnase-free  
Non-pyrogenic  
Filter and non-filter for choice  
Bag and rack are available to pack

## SPECIFICATIONS

Model	BPIC-707
Volume	200 µl
Unit	20



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)