



## STACKABLE 300µL TIP BOX BPIC-708

## STACKABLE 300?L TIP BOX BPIC-708

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## BPIC-708 STACKABLE 300?L TIP BOX



- Made from medical grade Polypropylene
- Slippy surface with excellent design
- DNase-free
- Rnase-free
- Non-pyrogenic
- Filter and non-filter for choice
- Bag and rack are available to pack

## SPECIFICATIONS

|        |             |
|--------|-------------|
| Model  | BPIC-708    |
| Volume | 300 $\mu$ l |
| Unit   | 20          |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)