



# 1000?L TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-715

## 1000µL TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-715

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

### BPIC-715 1000µL TIP BOX WITH ULTRA LOW RETENTION TIP



Made from medical grade Polypropylene  
Slippy surface with excellent design  
DNase-free  
Rnase-free  
Non-pyrogenic  
Filter and non-filter for choice  
Bag and rack are available to pack

## SPECIFICATIONS

Model	BPIC-715
Volume	1000 µl
Unit	10



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)