



# 5000?L TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-716

## 5000µL TIP BOX WITH ULTRA LOW RETENTION TIP BPIC-716

It offers quick and easy access to pipette tips, saving valuable time. Space-saving and environmentally compatible. Commonly used in liquid handling and molecular biology.

Used in Institutes, Pharmaceutical, Industrial, Microbiology, Plasmid Purification, Chemistry, Biology, Medical, Laboratory, Research, Liquid Handling Task.

## BPIC-716 5000µL TIP BOX WITH ULTRA LOW RETENTION TIP



- Made from medical grade Polypropylene
- Slippy surface with excellent design
- DNase-free
- Rnase-free
- Non-pyrogenic
- Filter and non-filter for choice
- Bag and rack are available to pack

## SPECIFICATIONS

Model	BPIC-716
Volume	5000 µl
Unit	10



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)