



## MANUAL PIPETTE CONTROLLER BPIP-407

## MANUAL PIPETTE CONTROLLER BPIP-407

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in laboratory and for research purpose.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Pipetaid, Pipet Gun, Laboratory Pipette Controller.

## BPIP-407 MANUAL PIPETTE CONTROLLER



- Chemical resistance
- Human engineering design
- Lightweight and Structure strong
- Lightweight and dexterous for operation
- A wide range volume for liquid controlling

## SPECIFICATIONS

Model	BPIP-407
Volume	0.1ml-100 ml
Colour	Yellow parts + Black cylinder
Unit	20



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)