



## MANUAL PIPETTE CONTROLLER BPIP-410

## MANUAL PIPETTE CONTROLLER BPIP-410

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in laboratory and for research purpose.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Pipetaid, Pipet Gun, Laboratory Pipette Controller.

## BPIP-410 MANUAL PIPETTE CONTROLLER



- Chemical resistance
- Human engineering design
- Lightweight and Structure strong
- Lightweight and dexterous for operation
- A wide range volume for liquid controlling

## SPECIFICATIONS

Model	BPIP-410
Volume	0.1ml-100 ml
Colour	Gray parts + Black cylinder
Unit	20



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)