



## MOTORIZED PIPETTE CONTROLLER BPIP-411

## MOTORIZED PIPETTE CONTROLLER BPIP-411

Pipette Controller offers an unmatched control of pipetting speeds from drop by drop or gravity dispense to fast liquid displacement. Designed for comfortable and uncomplicated pipetting even during long lasting sessions. It is light weight and easy to use in laboratory and for research purpose.

Used in Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology, Chemistry, Biology, Medical, Laboratory, Research.

Also known as Pipetaid, Pipet Gun, Laboratory Pipette Controller.

## BPIP-411 MOTORIZED PIPETTE CONTROLLER



- Chemical resistance
- Human engineering design
- Lightweight and Structure strong
- Lightweight and dexterous for operation
- A wide range volume for liquid controlling

## SPECIFICATIONS

Model	BPIP-411
Volume	0.1-100 ml
Colour	Dark blue

## OPTIONAL ACCESSORIES

Accessory Code	Name	MOQ
5300916006	RECHARGING STAND WALL MOUNT BRACKET	20
5300916007	CHARGER American plug	20
5300916008	CHARGER Europe plug	20
5300916009	CHARGER Australian plug	20
5300916010	EACH MEMBRANE FILTERS 0.45UM	20
5300916011	Battery	20



**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)