



## THERMAL STABILITY TESTER BPTL-209

## THERMAL STABILITY TESTER BPTL-209

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

## BPTL-209 THERMAL STABILITY TESTER



Eco-friendly metal bath heating, high effect of heat insulation. Air filling and releasing device with pressure gauge, adjusted by needle valve.

## SPECIFICATIONS

Model	BPTL-209
Working units	1 - 4
Temperature	Upto 400 °C
Metal bath heating	Eco-friendly, high effect of heat insulation
Test bomb	Standard stainless steel with inlet valve of 20m1
Safety protection device	Air filling and releasing device with pressure gauge, adjusted by needle valve.



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)