



CORROSION TESTER BPTL-224

CORROSION TESTER BPTL-224

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

BPTL-224 CORROSION TESTER



Sputtering multi-channel liquid circulator can guarantee temperature uniformity very well. Condenser quick detachable device

SPECIFICATIONS

| Model | BPTL-224 |
|--------------------------|---|
| Working units | 3, can be extended to 4-6 units |
| Circulator | Sputtering multi-channel liquid circulator can guarantee temperature uniformity |
| Condenser | Quick detachable device |
| Safety protection device | Can prevent over-low liquid level |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com