



**Product Image Coming Soon**

## **DISTILLATION TESTER BPTL-245**

## DISTILLATION TESTER BPTL-245

Petroleum testing is the analysis during upstream, midstream, and downstream production processes of petroleum products. It is most commonly used to test petroleum product, its product components, byproducts of crude oil, fuel, natural gas, upstream oil and gas and other formats of petroleum.

Used in Petroleum Industry, PVC Pipe Industry.

Also known as Distillation Apparatus.

## BPTL-245 DISTILLATION TESTER



Product Image Coming Soon

The instrument is suitable to determine the distillation characteristics of gasoline, aviation gasoline, jet fuels, and solvent having special boiling point, naphtha, diesel oil, distillate fuels and similar petroleum products. It adopts a special heating furnace to ensure the safety. The heating power can be adjusted continuously. The rising up and lowering down of distillation flask can be adjusted free by lifting device. Flexible and convenient to operate.

## SPECIFICATIONS

Model	BPTL-245
Distillation flask	125 ml
Water bath temperature controlling	
Range	(ambient temp.+10) °C ~60°C
Accuracy	±0.5°C
Receiving cylinder	100ml, scale division 1ml
Ambient temperature	15°C ~35°C
Relative humidity	≤ 85%
Power supply	AC(220±10%)V, 50Hz
Total power consumption	≤2000W
Dimension	460×400×550 mm
Net weight	18.5kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)