



## PURE WATER SUPPLY SYSTEM BPWS-102

## PURE WATER SUPPLY SYSTEM BPWS-102

This system can be combined with any brand's water purification system to supply pure water for specific equipment and it is easy to install and maintain. It has Independent power control and it can be run automatically according to the liquid level of the tank.

Used in Laboratory, Manufacturing, Reefkeeping, Aquarium, Laboratory, Research.

Also known as Laboratory Pure Water Supply.

## BPWS-102 PURE WATER SUPPLY SYSTEM



Independent power control and automatic operation, easy to install, use and maintain. Integrate level control, pressure pump, buffer tank and inlet valves together. It is unnecessary to connect to the circuit of pure water main-body. It can run automatically according to liquid level of the tank. Optional UV lamp module, to restrain bacteria's increase and reduce TOC.

## SPECIFICATIONS

Model	BPWS-102
Shape & Volume of tank	Roundness: 500, 1000, 2000 L Square: 200 L
Specifications of Pump	
Volume of Buffer Tank	1 L
Vertical self-priming height	2 meters
The Rated Stop Pressure	0.3 Mpa/0.4 Mpa, 0.35 Mpa/0.4 Mpa
Max Pressure/Pump Lift	About 0.3 Mpa/30 meters, About 0.35 Mpa/35 meters



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)