



## BENCHTOP SHAKING INCUBATOR BSBT-2201

# BENCHTOP SHAKING INCUBATOR BSBT-2201

## BSBT-2201 BENCHTOP SHAKING INCUBATOR



- Large LCD screen to display more data at same time
- R134a refrigerant, imported compressor and fan motor
- Big observation windows
- 304 Stainless steel chamber and platform, easy to clean
- The parameters can be automatically stored in case of power failure, and it will continue run as presetting program after turn on
- Microprocessor PID controller for temperature and shaking speed with timing function

## SPECIFICATIONS

Model	BSBT-2201
Speed Range	40-300 rpm
Electrical Requirement	220 V 50 Hz
Amplitude	26 mm
Display Resolution	0.1°C
Power Consumption	1050 W
Platform Dimensions	750x460 mm
External Dimension	1080Wx620Dx915H mm
Timing range	0~5999 min



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)