



BENCHTOP SHAKING INCUBATOR BSBT-2203

BENCHTOP SHAKING INCUBATOR BSBT-2203

BSBT-2203 BENCHTOP SHAKING INCUBATOR



- Large LCD screen to display more data at same time
- R134a refrigerant, imported compressor and fan motor
- Big observation windows
- 304 Stainless steel chamber and platform, easy to clean
- The parameters can be automatically stored in case of power failure, and it will continue run as presetting program after turn on
- Microprocessor PID controller for temperature and shaking speed with timing function

SPECIFICATIONS

Model	BSBT-2203
Speed Range	40-300 rpm
Electrical Requirement	220 V 50 Hz
Amplitude	26 mm
Display Resolution	0.1°C
Power Consumption	1300 W
Platform Dimensions	750x460 mm
External Dimension	1080Wx620Dx915H mm
Timing range	0~5999 min



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com