



BENCHTOP SHAKING INCUBATOR BSBT-301

BENCHTOP SHAKING INCUBATOR BSBT-301

Bench top shaking incubator has a compact design that saves your precious lab space. Stainless steel interiors with exceptional durability and trouble-free cleaning require less maintenance. Large transparent lid provides easy observation. Microprocessor controller permits easy operation. A variety of platforms and other lab accessories allows configuration to meet your demands.

Used in Cell Cultures, Cell Aeration, Microbiology, Increasing Solubility Rates, Metabolism Studies, Bacterial Cultures, Bacteriology.

Also known as Laboratory Incubator Shaker, Incubator Shaker.

BSBT-301 BENCHTOP SHAKING INCUBATOR



Speed can be set as 60-400 rpm

The electronic control of closed cycle adjustable heating system

Sainless steel interior and Electrostatic sprayed exterior

Call to restore function, without power supply

With strong influence of intermittent rapid heating system of.200 type and 100B type with air pressure regulating valve

Flexible and convenient, wide variety of detachable accessories

SPECIFICATIONS

Model	BSBT-301
Speed Range	20-300 rpm
Standard Capacity	500 mlx3, 250 mlx3, 100 mlx4, 50 mlx4
Temperature	RT+5°C~60°C
Temperature Accuracy	±0.1°C
Temperature Uniformity	±1°C
Shaking Speed Accuracy	±1 rpm
Timing Range	Optional
Motion	Orbital
Amplitude	24 mm
Tray	1
Platform Dimension	4375x405 mm
Inner Dimension	428x456x360 mm
Overall Dimension	640x605x520 mm
Display	LCD
Weight	72 kg
Power	365 W
Power Supply	220 V±10%



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com