



## FLOOR STANDING SHAKING INCUBATOR BSFT-2301

# FLOOR STANDING SHAKING INCUBATOR BSFT-2301

## BSFT-2301 FLOOR STANDING SHAKING INCUBATOR



Large LCD screen to display more data at same time  
R134a refrigerant, imported compressor and fan motor  
Big observation windows  
304 Stainless steel chamber and platform, easy to clean.  
There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement  
The parameters can be automatically stored in case of power failure, and it will continue run as presetting program after turn on  
Microprocessor PID controller for temperature and shaking speed with timing function

## SPECIFICATIONS

Model	BSFT-2301
Speed Range	40-300 rpm
Electrical Requirement	220 V 50 H
Amplitude	26 mm
Temperature Range	RT+5 - 65°C
Display Resolution	0.1°C
Timing Range	0~5999 min
Power Consumption	1900 W
Platform Dimension	750x460 mm
Platform quantity	2 pcs



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)