



## FLOOR STANDING SHAKING INCUBATOR BSFT-2303

# FLOOR STANDING SHAKING INCUBATOR BSFT-2303

## BSFT-2303 FLOOR STANDING SHAKING INCUBATOR



- Large LCD screen to display more data at same time
- R134a refrigerant, imported compressor and fan motor
- Big observation windows
- 304 Stainless steel chamber and platform, easy to clean.
- There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement
- The parameters can be automatically stored in case of power failure, and it will continue run as presetting program after turn on
- Microprocessor PID controller for temperature and shaking speed with timing function

## SPECIFICATIONS

|                        |            |
|------------------------|------------|
| Model                  | BSFT-2303  |
| Speed Range            | 40-300 rpm |
| Electrical Requirement | 220 V 50 H |
| Amplitude              | 26 mm      |
| Temperature Range      | 4 - 65°C   |
| Display Resolution     | 0.1°C      |
| Timing Range           | 0~5999 min |
| Power Consumption      | 2250 W     |
| Platform Dimension     | 750x460 mm |
| Platform quantity      | 2 pcs      |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)