



FLOOR STANDING SHAKING INCUBATOR BSFT-2304

FLOOR STANDING SHAKING INCUBATOR BSFT-2304

BSFT-2304 FLOOR STANDING SHAKING INCUBATOR



Large LCD screen to display more data at same time
R134a refrigerant, imported compressor and fan motor
Big observation windows
304 Stainless steel chamber and platform, easy to clean.
There is a 25mm instruction connection hole on the left side of the chamber for easy testing operation and temperature measurement
The parameters can be automatically stored in case of power failure, and it will continue run as presetting program after turn on
Microprocessor PID controller for temperature and shaking speed with timing function

SPECIFICATIONS

Model	BSFT-2304
Speed Range	40-300 rpm
Electrical Requirement	220 V 50 H
Amplitude	26 mm
Temperature Range	4 - 65°C
Display Resolution	0.1°C
Timing Range	0~5999 min
Power Consumption	2250 W
Platform Dimension	920x500 mm
Platform quantity	2 pcs



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com