



RECIPROCAL SHAKER BSRP-103

RECIPROCAL SHAKER BSRP-103

Light weight and portable reciprocal shaker with mild to vigorous shaking suitable for a variety of applications. Integral electronic timer handles both timed and continuous operations. Low noise and stable operation with over speed protection make it an ideal choice in different labs. Interchangeable platforms enable the use of different vessels. Used in Industry, Cosmetic, Pharmaceutical, Electronics, Laboratory.

Also known as Laboratory Shaker.

BSRP-103 RECIPROCAL SHAKER

Compact space-saving design fits easily in small spaces
Brushless DC motor ensures maintenance-free, reliable and quiet operation
Electronic time switching clock controls time or continuous operation mode
Simultaneous display of speed, timer and operating mode
RS232 interface enables data logging & computer control
Over speed detection and protection
Accommodate a variety of platforms, non-slip rubber mat, conical flasks and other sample containers

SPECIFICATIONS

Model	BSRP-103
Speed Range	100-500 rpm (orbital) 100-350 rpm (Linear)
Permissible Ambient Temperature	5-40°C
Shaking Orbit	10 mm
Permissible Relative Humidity	80.00%
Timing Range	1-1199 min
Motion	Orbital/Linear
Motor rating Input [W]	28
Motor rating output [W]	15
Maximum Load Capacity	7.5 kg
Operation	Continuous/Timed
Overall Dimension	420x370x100 mm
Display	LCD
Weight	13.5 kg
Power	30 W
Power Supply	110 V-240 V, 50/60 Hz



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com