



## THERMO SHAKER INCUBATOR BSTH-205

# THERMO SHAKER INCUBATOR BSTH-205

Used in Incubation, Cultivation, Re-suspending Pellets, Denaturin, PCR preparations, ELISA Assays, Colorimetric Assays, Enzyme, Protein Analysis.

Also known as Thermo Shaker, Lab Shaker, Thermo Shaker Incubator, Laboratory Thermo Shaker Incubator.

## BSTH-205 THERMO SHAKER INCUBATOR



Optimized 2D mixing designed for simultaneous mixing and vortexing applications

Pre-programmed soft keys and reliable universal holder

Automatic imbalance detection and anti-spill technology

Mixing in a 3-in-1 format: Plates (up to 96-well and 384-well formats), Micro-Tubes, Vortexing numerous tube formats

Custom blocks available for different experiment needs

## SPECIFICATIONS

Model	BSTH-205
Speed Range	300~2000 rpm
Temperature	RT+5-100°C
Temperature Accuracy	0.5°C
Temperature Resolution	0.1°C
Vortex Shaking Speed	None
Speed Increment	10 rpm
Shaking Orbit	3 mm
Time of Heating	≤12 min
Time of Cooling	No cooling function
Timing Range	1min ~ 99h59 min
Motion	Orbital
Overall Dimension	230x170x165 mm
Weight	5.5 kg
Power	150 W
Power Supply	230 V, 60 Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)