



## BIOSAFETY TRANSPORT BOX BTBT-105

## BIOSAFETY TRANSPORT BOX BTBT-105

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc.

Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

## BTBT-105 BIOSAFETY TRANSPORT BOX



Real time temperature value display

Refrigeration and temperature control to ensure the consistency of samples

Easy to operate and could be used repeatedly

Meet the basic requirements of cold chain transportation and storage of vaccines

## SPECIFICATIONS

Model	BTBT-105
Inner Capacity	20 L
Outer size	450x330x300 mm
Product Configuration	1 box 1 strap 1 digital thermometer 4 95KPa biosafety transport tanks 4 sponge reagent tubes 4 pieces of absorbent cotton 9 crystal ice boxes 1 partition 1 tray
Material	PP
Insulation Material	EPS/PU
Gross weight	6 kg
Packing size	490x350x330 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)