



BIOSAFETY TRANSPORT BOX BTBT-201

BIOSAFETY TRANSPORT BOX BTBT-201

The Biosafety Transport Box is a non-electric cooler used for the transportation of biological products requiring physical insulation, such as biological samples, pathogenic microorganism (virus) species, blood, etc.

Used in It can be used for the transportation of biological products requiring physical insulation, such as un2814, un2900, un3373 biological samples, pathogenic microorganism (virus) species, blood, etc.

BTBT-201 BIOSAFETY TRANSPORT BOX



AC & DC available Able to used in cars

SPECIFICATIONS

Model	BTBT-201
Capacity	30 L
Temperature Range	-25 to 30°C
Alarm	High&low temp, control error, power failure
External Dimension	610x360x460 mm
Power	DC12V; AC100~220V, 50/60Hz
Optional Parts	Basket; Datalogger



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com