



## 3D ROTATOR BTUR-301

## 3D ROTATOR BTUR-301

The versatile, end-over-end tube roller and Rotators are designed for mixing samples in both a vertical and horizontal plane. They can be placed in an environment temperature ranges from 5 to 40° and the rotation speed is adjustable, from 0 to 99 rpm. Tube Rollers, Rotators and Rocker are designed to deliver salient performance, accuracy, and reproducible results for all your application needs.

Used in Industrial, Clinical Laboratories, Academic Facilities, Environmental Applications, Food Technology, Wastewater Facilities.

Also known as Tube Rotator, Rotator, Laboratory Rotator.

## BTUR-301 3D ROTATOR



Unique combination of adjustable rotating frame.

A variety of rotation modes, from mild type to somersault type rotation mode.

Positive and negative rotation just gently press the rotating frame when instrument is working.

Widely used in molecular biology, histo-chemistry, bio-chemistry, clinical applications, etc.

Provide three types of rotating frame.

## SPECIFICATIONS

Model	BTUR-301
Speed Range	18 rpm
Voltage	AC 220V/AC 110V, 50/60Hz
Power	4 W
Dimension	240x145x210 mm
Net Weight	1 kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)