



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-701

# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-701

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-701 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



- Standard basket and lid
- Time and temperature real time display
- Fully microprocessor controller and programmable
- Digital display, easy operation for time and temperature
- Low level water protection
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

## SPECIFICATIONS

Model	BULC-701
Capacity	1.3 L
Ultrasonic frequency	40 KHz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	60 W
Heating power	200 W
Tank size	150x140x65 mm
Packing Size	260x260x260 mm
Gross Weight	2.5 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)