



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-704

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-704

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-704 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



Standard basket and lid

Time and temperature real time display

Fully microprocessor controller and programmable

Digital display, easy operation for time and temperature

Low level water protection

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

| | |
|--------------------------|----------------------|
| Model | BULC-704 |
| Capacity | 6 L |
| Ultrasonic frequency | 40 Khz |
| Time range/ Heating temp | 1-30 mins / 20- 80°C |
| Ultrasonic power | 180 W |
| Heating power | 400 W |
| Tank size | 300x150x150 mm |
| Packing Size | 430x280x380 mm |
| Gross Weight | 6 Kg |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com