



SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-708

SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-708

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-708 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



- Standard basket and lid
- Time and temperature real time display
- Fully microprocessor controller and programmable
- Digital display, easy operation for time and temperature
- Low level water protection
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

Model	BULC-708
Capacity	19 L
Ultrasonic frequency	40 KHz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	420 W
Heating power	600 W
Tank size	330x300x200 mm
Packing Size	460x430x430 mm
Gross Weight	11 Kg



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada
Email: contact@biolabscientific.com | Website: www.biolabscientific.com