



# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-710

# SINGLE FREQUENCY TYPE ULTRASONIC CLEANER BULC-710

Single Frequency type Ultrasonic Cleaner uses frequency of 40KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-710 SINGLE FREQUENCY TYPE ULTRASONIC CLEANER



- Standard basket and lid
- Time and temperature real time display
- Fully microprocessor controller and programmable
- Digital display, easy operation for time and temperature
- Low level water protection
- Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

## SPECIFICATIONS

Model	BULC-710
Capacity	30 L
Ultrasonic frequency	40 KHz
Time range/ Heating temp	1-30 mins / 20- 80°C
Ultrasonic power	600 W
Heating power	1000 W
Tank size	500x300x200 mm
Packing Size	630x430x430 mm
Gross Weight	15 Kg



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)