



DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER BULC-815

DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER BULC-815

Double adjustable frequencies ultrasonic cleaner uses frequency of 45 KHz /80 KHz with strong penetrating power and an appropriate cleaning solvent to clean items.

Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-815 DOUBLE ADJUSTABLE FREQUENCIES ULTRASONIC CLEANER

Over- Voltage and over current protection

Standard basket and lid

Low level water protection

Fully microprocessor controller and programmable

Time and temperature real time display

Standard power supply AC220+10%, 50/60Hz: 110V+10%, 60Hz

SPECIFICATIONS

| Model | BULC-815 |
|---------------------------------|----------------------|
| Capacity | 8 L |
| Adjustable ultrasonic frequency | 45 KHz /80 KHz |
| Time range/ Heating temp | 1-480 mins / RT-80°C |
| Ultrasonic power | 200 W |
| Heating power | 400 W |
| Tank size | 300x150x180 mm |
| Packing Size | 500x350x480 mm |
| Gross Weight | 15.5 Kg |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com