



## ULTRASONIC CLEANER BULC-901

## ULTRASONIC CLEANER BULC-901

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-901 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe  
The shell and inner cavity are made of 304 stainless steel, which is anti corrosion  
LED display

## SPECIFICATIONS

Model	BULC-901
Capacity	6.5 L
Ultrasonic frequency	Single frequency, 40 KHz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	180 W
Heating power	500 W
Tank size	300x150x150 mm
External Size	330x180x310 mm
Packing Size	430x250x350 mm
Gross weight	6 Kg
Power Supply	AC 220V, 50Hz



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)