



## ULTRASONIC CLEANER BULC-903

## ULTRASONIC CLEANER BULC-903

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-903 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe  
The shell and inner cavity are made of 304 stainless steel, which is anti corrosion  
LED display

## SPECIFICATIONS

|                          |                          |
|--------------------------|--------------------------|
| Model                    | BULC-903                 |
| Capacity                 | 15 L                     |
| Ultrasonic frequency     | Single frequency, 40 KHz |
| Time range/ Heating temp | 1-99 min/RT-80°C         |
| Ultrasonic power         | 360 W                    |
| Heating power            | 500 W                    |
| Tank size                | 330x300x150 mm           |
| External Size            | 360x330x310 mm           |
| Packing Size             | 460x400x400 mm           |
| Gross weight             | 10 Kg                    |
| Power Supply             | AC 220V, 50Hz            |



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada  
Email: [contact@biolabscientific.com](mailto:contact@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)