



ULTRASONIC CLEANER BULC-910

ULTRASONIC CLEANER BULC-910

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

BULC-910 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion
LED display

SPECIFICATIONS

| | |
|--------------------------|--------------------------|
| Model | BULC-910 |
| Capacity | 30 L |
| Ultrasonic frequency | Single frequency, 40 KHz |
| Time range/ Heating temp | 1-60 min/RT-80°C |
| Ultrasonic power | 600 W |
| Heating power | 500 W |
| Tank size | 300x150x150 mm |
| External Size | 530x330x350 mm |
| Packing Size | 640x400x400 mm |
| Gross weight | 15 Kg |
| Power Supply | AC 220V, 50Hz |



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com