



## ULTRASONIC CLEANER BULC-914

## ULTRASONIC CLEANER BULC-914

Ultrasonic Cleaner is a multi purpose cleaning device which uses mechanical vibrations for appropriate cleaning solvents Used in Pharmaceutical, Chemical and physical cleaning, electronics industry, Research, laboratory. Also known as Digital ultrasonic cleaner, Ultrasonic cleaning bath, Ultrasonic cleaning tank.

## BULC-914 ULTRASONIC CLEANER



Standard accessories: soundproof cover, basket, PVC pipe

The shell and inner cavity are made of 304 stainless steel, which is anti corrosion  
LED display

## SPECIFICATIONS

Model	BULC-914
Capacity	22.5 L
Ultrasonic frequency	Double frequency 28/40 Khz
Time range/ Heating temp	1-99 min/RT-80°C
Ultrasonic power	480 W
Heating power	500 W
Tank size	500x300x150 mm
External Size	530x330x310 mm



**Biolab Scientific Ltd.**

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)