



MINI MIXER BVMX-303

MINI MIXER BVMX-303

This sturdy, reliable mixers can adjust speed to meet different experimental requests. Multiple mixing modes, Touch and continuous operations. Aluminium cast base avoids unnecessary "walking". A variety of optional adapters are available which process 70 samples at the same time to make the experiment more convenient and efficient. Wide speed range:0-3000rpm. Built-in short mix, Microprocessor control, LED display the speed and time.

Used in Mixing, Laboratory, Research, Medical.

Also known as Laboratory Vortex Mixer.

BVMX-303 MINI MIXER



Mini and compact appearance, powerful function, sucker type under-chassis, super shockproof, ideal for high speed operation

Multiple mixing modes, touch-operation, continuous operation and timing set

Speed adjustable, and wide speed range, Stable low speed, powerful high speed

Shaking blocks convenient for replacement, stable and reliable, sturdy and durable eccentric bearing designed

Various tube holders are for optional (Suitable for Eppendorf Tubes)

SPECIFICATIONS

Model	BVMX-303
Speed Range	300 ~ 3000 rpm
Shaking Orbit	3 mm
Shaking Mode	Circle
Operation	key + digital Display
Operation Mode	continuous running or short mix
Voltage	DC 24 V
Dimension	W.132 x D.156 x H.150 mm
Net Weight	2.2 kgs
Power	30 W

ACCESSORIES

Accessory Code	Name	Description
3100808021	Standard head	Suit for tubes and containers with dia.<30 mm
3100808022	Standard platform	Standard platform to hold foam tube holder
3100808023	Convex pad	Use with MX-P. suit for tubes and containers with dia.<99 mm

OPTIONAL ACCESSORIES

Accessory Code	Name	DESCRIPTION
3100808006	Standard head	Suit for tubes and containers with dia.<30 mm

Accessory Code	Name	DESCRIPTION
3100808014	Tube holder platform	for 1.5/2.0 ml tubes x 22

3100808007	Standard platform Foam tube holder	Standard platform to hold foam tube holder
3100808008	Convex pad	Use with MX-P. suit for tubes and containers with dia.<99 mm
3100808009	Foam tube holder	for 0.2 ml tubes 70 holes. use with MX-P
3100808010	Foam tube holder	for 0.5 ml tubes 36 holes. use with MX-P
3100808011	Foam tube holder	for 1.5/2.0 ml tubes 20 holes. use with MX-P
3100808012	Foam tube holder	for 15 ml tubes 8 holes. use with MX-P
3100808013	Foam tube holder	for 50 ml tubes 4 holes. use with MX-P

3100808015	Tube holder platform	for 15 ml tubes x 4
3100808016	Tube holder platform	for 50 ml tubes x 2
3100808017	Plate holder platform for all standard plates	for all standard plates
3100808018	Universal platform	for tube stands below and all standard plates
3100808019	Form tube holder fixer	for MX-C & MX-D & MX-E
3100806020	Form tube holder fixer	for MX-F & MX-G



Biolab Scientific Ltd.

3660 Midland Avenue, Suite 300, Toronto, Ontario M1V 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com