



HOT AIR STERILIZER

HOT AIR STERILIZER BFO1B1 BFO1B2

PORTABLE PRESSURE STEAM STERILIZER

Durable, corrosion resistant and easy to clean stainless steel design in a wide range of capacities. Our hot air sterilizers are microprocessor controlled and perform reliably. A digital display with a broad observation window allows easy operation.



SPECIFICATIONS

Model	BFO1B1	BFO1B2
Old Model	BSHA-207	BSHA-208
Capacity	18 L	24 L
Function	Electric or LPG heated.	
Sterilization Temperature	126°C	
Working Pressure	0.14-0.16 Mpa	
Maximum Working Pressure	0.165 Mpa	
Chamber Dimension	280x260 mm	280x390mm
Overall Dimension (LxWxH)	410x410x430 mm	410x410x550 mm
Chamber Material	Stainless Steel	
Weight	14/16 kg	15/17 kg
Power	2 KW	
Power Supply	220 V, 50 Hz	
Alt Name	Portable Pressure Steam Sterilizer	

APPLICATIONS

Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research.

HOT AIR STERILIZER BFO1C1 BFO1C2

PORTABLE PRESSURE STEAM STERILIZER

Durable, corrosion resistant and easy to clean stainless steel design in a wide range of capacities. Our hot air sterilizers are microprocessor controlled and perform reliably. A digital display with a broad observation window allows easy operation.



SPECIFICATIONS

Model	BFO1C1	BFO1C2
Old Model	BSHA-209	BSHA-210
Capacity	18 L	24 L
Function	Electric heated	
Sterilization Temperature	126°C	
Steam Exhausting Temperature	≤1°C	
Working Pressure	0.14-0.16Mpa	
Maximum Working Pressure	0.165 Mpa	
Chamber Dimension	280x260 mm	280x390mm
Overall Dimension (LxWxH)	410x410x520 mm	410x410x630 mm
Chamber Material	Stainless Steel	
Weight	14/16 kg	15/17 kg
Power	2 KW	
Power Supply	220 V, 50 Hz	
Alt Name	Portable Pressure Steam Sterilizer	

APPLICATIONS

Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research.

HOT AIR STERILIZER BFO1D1 BFO1D2

PORTABLE PRESSURE STEAM STERILIZER

Durable, corrosion resistant and easy to clean stainless steel design in a wide range of capacities. Our hot air sterilizers are microprocessor controlled and perform reliably. A digital display with a broad observation window allows easy operation.

SPECIFICATIONS

Model	BF01D1	BF01D2
Old Model	BSHA-211	BSHA-212
Capacity	18 L	24 L
Function	Electric heated	
Sterilization Temperature	126°C	
Steam Exhausting Temperature	≤1°C	
Working Pressure	0.14-0.16Mpa	
Maximum Working Pressure	0.165 Mpa	
Timer Range	0-99min	
Chamber Dimension	280x260mm	270x500mm
Overall Dimension (LxWxH)	440x420x630 mm	440x420x750 mm
Chamber Material	Stainless Steel	
Weight	12/14 kg	14/15 kg
Power	2 KW	
Power Supply	220 V, 50 Hz	
Alt Name	Portable Pressure Steam Sterilizer	

APPLICATIONS

Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com