



CUBE ICE MAKER

CUBE ICE MAKER BFC1D1 TO BFC1D4

SPLIT-TYPE ICE MAKER



Split-type structure, large ice bin, large ice storage capacity and strong applicability;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Square ice cubes with beautiful shapes.
Square

SPECIFICATIONS

Model	BFC1D1	BFC1D2	BFC1D3	BFC1D4
Old Model	BIFL-213	BIFL-214	BIFL-215	BIFL-216
Ice Making Capacity	200kg/24h	300kg/24h	500kg/24h	1000kg/24h
Ice Storage Capacity	100kg	150kg		500kg
Ice Shape	Square			
Power/Voltage (W/V)	990w/220v	1240w/220v	2300w/220v	4500w/220v
Cooling Type	Water cooling			
Product Dimension (mm)	560x800x1530	760x850x1620	760x850x1890	1250x950x2000
Machine Appearance	Split machine			
Alt Name	Split-Type Ice Maker			



BFC1D1 BFC1D2 BFC1D3



BFC1D4

FEATURES

Split-type structure, large ice bin, large ice storage capacity and strong applicability;
Fully automatically computer controlled, convenient for use;
Automatic detection and protection in case of operation, water deficiency, ice full and malfunction;
Highly efficient fluoride-free compressor, energy efficient and environment-friendly;
Square ice cubes with beautiful shapes.



Square



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com