



PORTABLE REFRIGERATOR

PORTABLE REFRIGERATOR BEM2F2 TO BEM2F10

PORTABLE REFRIGERATOR

The Portable Refrigerator is a versatile cooling solution designed for travel, camping, medical, and laboratory use, providing precise and reliable temperature control in a compact, durable design. Ideal for storing vaccines, food, and other temperature-sensitive items, this refrigerator combines portability with energy efficiency to ensure stable cooling in any environment.



SPECIFICATIONS

Model	BEM2F2	BEM2F3	BEM2F4	BEM2F5
Old Model	BPOR-100	BPOR-101	BPOR-102	BPOR-103
Capacity	6L	8L	12L	15L
Internal&External Materials	PP+PE			
Insulation Material	PU			
Temp. Resolution	0.1°C			
Manufacturing Process	Blow molding+injection molding, burr-free surface			
Cooling Source Material	Sodium polyacrylate			
Cooling Retention Time	Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C			
Configuration	Refrigerator*1, phase-change ice box(400ml)*4, digital thermometer*1, shoulder strap*1	Refrigerator*1, phase-change ice box(400ml)*4, digital thermometer*1	Refrigerator*1, phase-change ice box(400ml)*6, digital thermometer*1, shoulder strap*1, tray*1	
Internal Size(W*D*H)	240*145*180mm	295*167*168mm	325*170*225mm	363*203*244mm
External Size(W*D*H)	330*235*275mm	335*225*205mm	425*262*305mm	
Package Size(W*D*H)	380*270*300mm		460*355*360mm	
Net Weight	2.25kg	2kg	3.5kg	3.7kg
Gross Weight	2.7kg	2.4kg	4kg	4.3kg
Alt Name	Portable Refrigerator			

Model	BEM2F6	BEM2F7	BEM2F8	BEM2F9	BEM2F10
Old Model	BPOR-104	BPOR-105	BPOR-106	BPOR-107	BPOR-108
Capacity	20L	33L	55L	65L	85L
Internal&External Materials	PP+PE				
Insulation Material	PU				
Temp. Resolution	0.1°C				
Manufacturing Process	Blow molding+injection molding, burr-free surface				
Cooling Source Material	Sodium polyacrylate				

Cooling Retention Time	Maintains an internal temperature of 2~8°C for over 24 hours at an ambient temperature up to 43°C				
Configuration	Refrigerator*1, phase-change ice box(400ml)*9, digital thermometer*1, shoulder strap*1, tray*1, divider panel*1	Refrigerator*1, phase-change ice box(400ml)*12, digital thermometer*1		Refrigerator*1, phase-change ice box(400ml)*15, digital thermometer*1	
Internal Size(W*D*H)	380*245*230mm	480*250*290mm	590*335*320mm	540*425*315mm	580*465*320mm
External Size(W*D*H)	450*330*300mm	570*310*360mm	640*395*395mm	620*485*375mm	640*520*375mm
Package Size(W*D*H)	490*360*340mm	585*325*390mm	670*430*430mm	640*505*395mm	650*540*400mm
Net Weight	3.8kg	4.5kg	6.3kg	7.3kg	8kg
Gross Weight	4.5kg	5.25kg	7.3kg	8kg	9.3kg
Alt Name	Portable Refrigerator				



BEM2F2



BEM2F3



BEM2F4 BEM2F5 BEM2F7



BEM2F6



BEM2F8



BEM2F9 BEM2F10



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com