



MULTISTATION STIRRER

MULTISTATION STIRRER BER1BR1 BER1BR2

FOUR ROWS MAGNETIC STIRRER



1. 4 unit stainless hotplate can work together or separately.
2. Stirring speed is stepless adjusted.
3. Magnetic mixing equipment.

SPECIFICATIONS

Model	BER1BR1	BER1BR2
Function		
Max. Volume (ml)	1000ml*4	3000ml*4
Speed (r.p.m)	100 ~ 2000	
Working Time	Continuous	
Minimum Resolution	100	
Structure		
Exterior	Cold plate spraying, stainless steel plate	
Driving Mode	Motor	
Motor Power	0.01kW	0.02kW
Controller		
Continuous Operation	Yes	
Setting Mode	Scale	
Speed Display Mode	Scale	
Timer	No	
Operation Function	No	
Specification		
Exterior Size (W*L*H)(mm)	300*270*80	400*320*120
Packing Size (W*L*H)(mm)	340*310*110	440*360*150
Plate Size (mm)	300*250	400*300
Plate Load	5KG	20KG
Power supply (50/60Hz) Current Rating	110 or 220V	
NW/GW (kg)	3/3.6	6/7
Alt Name	Four Rows Magnetic stirrer	



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com