



WATER DISTILLER

WATER DISTILLER BFG1E1

WATER DISTILLER



SPECIFICATIONS

Model	BFG1E1
Power	750W
Voltage	220V/50Hz
Capacity	1.5L/h
Chamber size	780mm*200mm
Package size	290mm*290mm*390mm
Product size	245mm*285mm*510mm
N.W.	3.5kgs
G.W.	5.0kgs
Alt Name	Water Distiller

WATER DISTILLER BER1CN1 BER1CN2 BER1CN3

STAINLESS STEEL WATER DISTILLER (AUTOMATIC CONTROL TYPE)



1. The whole shell is designed in stainless steel, and the evaporation tank, condenser and heater are all made of stainless steel, which are resistant to corrosion and aging.
2. With water shortage protection function, when the liquid level is too low, it will automatically alarm and quickly cut off the heating, to ensure that the heater will not burn out, it is easy to operate, safe and reliable.
3. Automatic water replenishment function, when the cooling water is restored to supply water, it will automatically enter the evaporation bucket, to ensure that the float switch automatically controls the power supply, and automatically restore the output of distilled water.

SPECIFICATIONS

Model	BER1CN1	BER1CN2	BER1CN3
Voltage(V)	220	380	
Output(L/H)	≥5	≥10	≥20
Heating Power(kw)	4.5	7.5	15

G.W(kg)	9.5	10.5	13.5
N.W(kg)	6	7.5	11
Exterior Size (DiaxH) mm	Φ 305x730	Φ 335x830	Φ 405x980
Packing Size (LxWxH) mm	380x350x870	390x350x900	410x410x1055
Alt Name	Stainless Steel Water Distiller		

APPLICATIONS

This product uses tap water as the source water and uses electric heating distillation method to produce distilled water, which is suitable for medical treatment, chemical industry, scientific research , laboratories, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com