



DIALYSIS BED

DIALYSIS BED BGX1E1

MANUAL 2 FUNCTION DIALYSIS BED

Dialysis Bed is mainly used for setting purposes during dialysis treatment. This bed is best to be used when patients with end-stage renal disease spend many hours each week undergoing hemodialysis. We contribute toward enhanced care, through our innovative range of seating solutions. Our specialised dialysis bed has been created with both patient comfort and nurse safety in mind.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock
Bed sheet top made from perforated steel, epoxy powder coated and baked finish, 3 sections Bed framework made from mild steel ERW rectangular tubes
With casters dia. 5" with individual brake
With equipment shelf on the foot end
2 functions adjusted by hidden type hand crank
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
Optional Accessories: Shoe shelf, mosquito net rack, dinner table.

SPECIFICATIONS

Model	BGX1E1
Old Model	BHBD-302
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500H mm
Alt Name	Manual 2 Function Dialysis Bed

APPLICATIONS

Dialysis treatment, Hemodialysis.

DIALYSIS BED BGX1F1

MANUAL 2 FUNCTION DIALYSIS BED

Dialysis Bed is mainly used for setting purposes during dialysis treatment. This bed is best to be used when patients with end-stage renal disease spend many hours each week undergoing hemodialysis. We contribute toward enhanced care, through our innovative range of seating solutions. Our specialised dialysis bed has been created with both patient comfort and nurse safety in mind.



ABS head & foot boards can be removable. with bumper wheels, 12kg/pair, color option: brown, blue, green, with patient name card holder
2 pcs Aluminium side rails with 6 S.S bars, up and down to protect the patient not fall down from the bed, collapsed type with safe lock
Bed sheet top made from perforated steel, epoxy powder coated and baked finish, 3 sections Bed framework made from mild steel ERW rectangular tubes
Central locking ABS foot pedal with casters dia. 5"
With equipment shelf on the foot end
2 functions adjusted by hidden type hand crank
Accessories: 1 pc SS telescopic I.V pole; 4 pcs urine bag hooks; 4 pcs I.V holes
Optional Accessories: Shoe shelf, mosquito net rack, dinner table.

SPECIFICATIONS

Model	BGX1F1
Old Model	BHBD-303
Function	2
Back-rest lifting angle	0°-80°±5°
Foot-rest lifting angle	0°-40°±5°
IV pole	Yes
Dimension	2160Lx980Wx500H mm
Alt Name	Manual 2 Function Dialysis Bed

APPLICATIONS

Dialysis treatment, Hemodialysis



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com