



## MANUAL HOSPITAL BED

## MANUAL HOSPITAL BED BGZ1V1

### FLAT HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



## SPECIFICATIONS

Model	BGZ1V1
Old Model	BHBD-521
Size	2150*900*500mm
Material	ABS plastic with punching
Alt Name	Flat Hospital Bed

## APPLICATIONS

Recuperating, Nursing

## MANUAL HOSPITAL BED BGZ1J1

### SINGLE CRANK HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical units and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



## SPECIFICATIONS

Model	BGZ1J1
Old Model	BHBD-522
Specification (L x W x H)	2150 x 960 x 500 mm
Back Adjustment	≥ 75°
Leg Adjustment	≥ 35°
Load Capacity	240 kg
Bed Frame	Made of 80 x 40 x 1.5 mm formed square tubes, robot-welded; strong structure capable of carrying ≥ 240 kg

Bed Surface	1.2 mm cold-rolled steel plate, stamped into concave shape with multiple pores for breathability and anti-slip function
Alt Name	Single Crank Hospital Bed

## APPLICATIONS

Recuperating, Nursing.

## MANUAL HOSPITAL BED BGZ1F1

CRANKS HOSPITAL BED



## SPECIFICATIONS

Model	BGZ1F1
Old Model	BHBD-523
Size	2150*960*500mm
Back section adjustment	≥75°
Loading capacity	≥240kg
Uraine hook	Yes
Bumpers	Yes
Alt Name	Cranks Hospital Bed


## ITEMS INCLUDED

No	Name	Quantity
1	Stainless steel handle	1 pc
2	Aluminum alloy Guardrails	2 pcs
3	Universal 125 castors with braking	4 pcs
4	Bedhead board	2 pcs
5	IV Pole	1 pc

## ACCESSORIES FOR PURCHASE

No	Name
1	Night light
2	Mattress
3	Overbed Food Table
4	Bedside cabinet
5	Color is optional

## FEATURES



## MANUAL HOSPITAL BED BGW1K1

MANUAL BED



## SPECIFICATIONS

Model	BGW1K1
Old Model	BHBD-509
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2080 mm (L) * 980 mm (W) * 470 mm (H)
Bed board dimension	1910 mm (L) * 830 mm (W)
Functions	Two functions controlled by cranks
Backrest	0 - 75°
Knee gatch	0 - 45°
Height (bed board to ground)	470 mm
Protective railings	Collapsible 5-column aluminum protective railings, easy to grasp
Head & foot boards	PP head & foot boards, integrally blow-molded, easy to detach and clean
Castors	360° swivel castors (Dia. 125 mm) with individual brakes
Optional feature	Shoe rack under bed board (optional)
Net weight	65 KG
Standard accessory	I.V. Pole
Alt Name	Manual Bed



## FEATURES



## MANUAL HOSPITAL BED BGW1L1

3 CRANKS MANUAL BED



## SPECIFICATIONS

Model	BGW1L1
Old Model	BHBD-510
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2160 mm (L) * 970 mm (W) * 575 mm - 750 mm
Bed board dimension	1930 mm (L) * 830 mm (W)
Functions	Three functions controlled by cranks
Backrest	0 - 65°
Knee gatch	0 - 45°
HI-LOW	750 mm - 575 mm
Protective railings	Collapsible 6-column aluminum protective railings, easy to grasp
Head & foot boards	PP head & foot boards, integrally blow-molded, easy to detach and clean
Castors	360° swivel castors (Dia. 125 mm) with individual brakes
Net weight	102 KG
Standard accessory	I.V. Pole
Alt Name	3 Cranks Manual bed

## FEATURES



## MANUAL HOSPITAL BED BGZ1L1

### MANUAL HOSPITAL BED WITH 2 FUNCTION

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



## SPECIFICATIONS

Model	BGZ1L1
Old Model	BHBD-525
Size (L x W x H)	2200 x 1020 x 420/780 mm

Back-rest Adjustment	65° (±10°)
Knee-rest Adjustment	35° (±10°)
Height Range	420-780 mm
Alt Name	Manual Hospital Bed With 2 Function

## ITEMS INCLUDED

No	Name	Quantity
1	European style guardrail	4 pcs
2	ABS head and foot board	2 pcs
3	Crank quantity	2 pcs
4	5" Castors with brake	4 pcs
5	Urine hook	1 pc
6	Protective bumpers	4 pcs
7	IV pole	1 pc

## APPLICATIONS

Recuperating, Nursing.

## MANUAL HOSPITAL BED BGZ1K1

### MANUAL HOSPITAL BED

Manual Hospital Bed is suitable for various operations conducted by medical unit and recommended when there is a caregiver available to make adjustments. As a cost-effective solution, these manual beds have the option to operate via use of hand cranks for adjusting the height of the bed, and the head and the feet as well. They also benefit from the lack of risk of motor failure and maintenance.



Material: ABS  
 4 pcs of European style guardrail, upside down structure  
 Luxury central control castors, brake the bed with one stamp  
 With 4 bumpers

## SPECIFICATIONS

Model	BGZ1K1
Old Model	BHBD-524
Size (L x W x H)	2280 x 1120 x 420/780 mm
Height Adjustment	420-780 mm
Back-rest Adjustment	65° (±10°)
Knee-rest Adjustment	35° (±10°)
Weight Loading Capacity	> 240 kg
Alt Name	Manual Hospital Bed

## ITEMS INCLUDED

No	Name	Quantity
1	European Style Guardrail	4 pcs
2	ABS Head & Foot Board	2 pcs

3	5" Castors with Brake	4 pcs
4	Urine Hook	1 pc
5	Protective Bumpers	4 pcs
6	IV Pole	1 pc
7	Crank Quantity	3 pcs

## APPLICATIONS

Recuperating, Nursing.

## MANUAL HOSPITAL BED BGW1M1

4 CRANKS MANUAL BED



## SPECIFICATIONS

Model	BGW1M1
Old Model	BHBD-511
Material	Powder coated steel bed frame, top class ABS & PP plastic
Overall dimension	2120 mm (L) * 1055 mm (W) * 505 mm - 665 mm
Bed board dimension	1930 mm (L) * 900 mm (W)
Functions	Five functions activated by cranks
Backrest	0 - 60°
Knee gatch	0 - 30°
Trend/Reverse Trend	0 - 7°
HI-LOW	665 mm - 505 mm
Head & foot boards	PP head & foot boards, integrally blow-molded, easy to detach and clean
Protective railings	4 damping PP protective railings, controlled by gas springs, seamless connection with stretcher
Castors	360° swivel castors (Dia. 150 mm) with central brake under foot board
Net weight	117 KG
Standard accessory	I.V. Pole
Alt Name	4 Cranks manual bed

## FEATURES







**Biolab Scientific Ltd.**

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada  
Email: [info@biolabscientific.com](mailto:info@biolabscientific.com) | Website: [www.biolabscientific.com](http://www.biolabscientific.com)