



ORTHOPEDICS TRACTOR RACK

ORTHOPEDICS TRACTOR RACK BGY1CE1

ORTHOPEDICS TRACTOR RACK

Biolab's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.



Adopt carbon fiber material, the perspective imaging is clear, light and strong
The plug in design is convenient and quick to connect with table
Equipped with a traction frame accessory car, which is convenient to move and save space

SPECIFICATIONS

| | |
|-----------------------------------|--------------------------|
| Model | BGY1CE1 |
| Old Model | BOPT-403 |
| Table Length | 1550±50 mm |
| Table Width | 450±50 mm |
| Traction device horizontal rotary | 0-180° |
| Traction distance | (0-200)+10% |
| Traction Device height adjustment | >190 mm |
| Alt Name | Orthopedics tractor rack |

APPLICATIONS

Surgery, Leg operation, Medical unit

ORTHOPEDICS TRACTOR RACK BGY1CC1

ORTHOPEDICS TRACTOR RACK

Biolab's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.



Adopt carbon fiber material, the perspective imaging is clear, light and strong
The plug in design is convenient and quick to connect with table
Equipped with a traction frame accessory car, which is convenient to move and save space

SPECIFICATIONS

| | |
|--------------|------------|
| Model | BGY1CC1 |
| Old Model | BOPT-401 |
| Table Length | 1420±50 mm |

| | |
|-----------------------------------|--------------------------|
| Table Width | 765±50 mm |
| Traction device horizontal rotary | 0-180° |
| Traction distance | (0-200)+10% |
| Traction Device height adjustment | >400 mm |
| Alt Name | Orthopedics tractor rack |

APPLICATIONS

Surgery, Leg operation, Medical unit

ORTHOPEDICS TRACTOR RACK BGY1CD1

ORTHOPEDICS TRACTOR RACK

Biolab's Orthopedics tractor rack is mainly applied to orthopedics surgical by using this tractor rack the patient's leg can be performed side traction, prone traction and open-up position of the body during the operation.



Adopt carbon fiber material, the perspective imaging is clear, light and strong
The plug in design is convenient and quick to connect with table
Equipped with a traction frame accessory car, which is convenient to move and save space

SPECIFICATIONS

| | |
|-----------------------------------|--------------------------|
| Model | BGY1CD1 |
| Old Model | BOPT-402 |
| Table Length | 1700±10 mm |
| Table Width | 750±100 mm |
| Traction device horizontal rotary | 0-180° |
| Traction distance | (0-200)+10% |
| Alt Name | Orthopedics tractor rack |

APPLICATIONS

Surgery, Leg operation, Medical unit



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com