



PIPETTE BULB

PIPETTE BULB BHS1B1

RUBBER PIPETTE BULB



Our pinch and release valve system allows for one-hand operation
Can accommodate pipettes up to 25mL
Easy to use
Made from natural rubber
Chemical resistant rubber filler with 3 valves for suction
Two valve system with a removable top valve for easy cleaning and one-handed use

SPECIFICATIONS

Model	BHS1B1
Volume	2-25 ml
Color	Red
Unit	50
Alt Name	Rubber Pipette Bulb

APPLICATIONS

Medical, Laboratory, Research, Chemistry, Biology, Liquid Handling Task, Institutes, Pharmaceutical, Industrial, Microbiology



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com