



STRONG ACID AND ALKALI STORAGE CABINET

STRONG ACID AND ALKALI STORAGE CABINET BEM2C1 TO BEM2C8

STRONG ACID & ALKALI STORAGE CABINET

The Strong Acid & Alkali Storage Cabinet is a specialized safety cabinet designed for the safe storage of corrosive chemicals, including strong acids and alkalis, in laboratory and industrial environments.



1. Padlock design.
2. Adjustable shelves.
3. Anti-corrosive materials for strong acid and alkali.
4. Customizable.

SPECIFICATIONS

Model	BEM2C1	BEM2C2	BEM2C3	BEM2C4
Old Model	BCBT-201	BCBT-202	BCBT-203	BCBT-204
Capacity	4 gallon(15L)	12 gallon(45L)	20 gallon(75L)	28 gallon(105L)
Material	PP			
Adjustable Shelf	Optional	1pc		2pcs
Door	Single/Manual		Double/Manual	
External Size(W*D*H)mm	430*430*560	590*460*890	910*590*890	910*600*900
Package Size(W*D*H)mm	515*515*695	675*545*1025	995*675*1025	995*685*1035
Gross Weight	16kg	30kg	46kg	48kg
Alt Name	Strong Acid & Alkali Storage Cabinet			

Model	BEM2C5	BEM2C6	BEM2C7	BEM2C8
Old Model	BCBT-205	BCBT-206	BCBT-207	BCBT-208
Capacity	30 gallon(114L)	45 gallon(170L)	60 gallon(227L)	90 gallon(340L)
Material	PP			
Adjustable Shelf	2pcs			
Door	Double/Manual			
External Size(W*D*H)mm	1090*460*1120	1090*460*1650	860*860*1650	1090*860*1650
Package Size(W*D*H)mm	1175*545*1255	1175*545*1785	945*945*1785	1173*943*1785
Gross Weight	58kg	72kg	90kg	107kg
Alt Name	Strong Acid & Alkali Storage Cabinet			



BEM2C1 BEM2C2



BEM2C3 BEM2C4 BEM2C5 BEM2C6



BEM2C7 BEM2C8



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada

Email: info@biolabscientific.com | Website: www.biolabscientific.com