



MECHANICAL ULTRASONIC DESKTOP CLEANER

MECHANICAL ULTRASONIC DESKTOP CLEANER BJR1L1 TO BJR1L10

MECHANICAL ULTRASONIC DESKTOP CLEANER



1. Mechanical temperature and time setting;
2. Industrial high Q transducers, cleaning conversion efficiency is higher;
3. No-stud welding technology, cleaning effect is better;
4. MCU-SWEEP ultrasonic generator frequency-drive, ultrasonic cleaning more uniform, more powerful;
5. MCH heating system, maximum temperature up to 80°C;
6. All stainless steel material, inner tank is made by punching SUS304 1.0mm;
7. 2500V high voltage test to ensure the product safety, strict production process, aging process and quality control
8. All products are CE, FCC, RoHS, PSE certificated.

SPECIFICATIONS

Model	BJR1L1	BJR1L2	BJR1L3	BJR1L4
Old Model	BULC-402	BULC-403	BULC-404	BULC-405
Capacity (L)	2.0	3.2	5.0	4.5
Frequency Working (KHz)	40KHz			
Ultrasonic power (W)	70	120		180
Heating power (W)	100		200	
Timing time (mins)	1~120min			
Temperature range (°C)	20-80			
Cover	YES			
Valve Drain	NO			YES
Washing Basket	optional			
Inner Tank Dimension (mm)	150x135x100	240x135x100	240x135x150	300x150x100
Packing Size (mm)	175x160x217	268x170x243	268x170x290	325x178x243
Gross weight (Kg)	2.7	4.77	5.4	5.65
Alt Name	Mechanical Ultrasonic Desktop Cleaner			

Model	BJR1L5	BJR1L6	BJR1L7	BJR1L8
Old Model	BULC-406	BULC-407	BULC-408	BULC-409
Capacity (L)	6.8	11	15	20
Frequency Working (KHz)	40KHz			
Ultrasonic power (W)	180	240	360	
Heating power (W)	200	300	400	
Timing time (mins)	1~120min			
Temperature range (°C)	20-80			
Cover	YES			
Valve Drain	YES			
Washing Basket	optional			
Inner Tank Dimension (mm)	302x150x150	300x240x150	330x300x150	330x300x200
Packing Size (mm)	325x178x290	368x270x290	397x320x317	397x320x368
Gross weight (Kg)	6.2	8.5	11.5	12.5

Alt Name	Mechanical Ultrasonic Desktop Cleaner	
Model	BJR1L9	BJR1L10
Old Model	BULC-410	BULC-411
Capacity (L)	22	30
Frequency Working (KHz)	40KHz	
Ultrasonic power (W)	480	600
Heating power (W)	500	
Timing time (mins)	1~120min	
Temperature range (°C)	20-80	
Cover	YES	
Valve Drain	YES	
Washing Basket	optional	
Inner Tank Dimension (mm)	500x300x150	500x300x200
Packing Size (mm)	560x320x317	560x320x385
Gross weight (Kg)	14.5	16.0
Alt Name	Mechanical Ultrasonic Desktop Cleaner	



BJR1L1



BJR1L2



BJR1L3



BJR1L4



BJR1L5 BJR1L6



BJR1L7



BJR1L8



BJR1L9



BJR1L10

APPLICATIONS

Ultrasonic cleaner are apply to the cleaning and disinfection of medical tools and instruments in the dental clinic. Others applications include the fast and efficient cleaning and disinfection/sterilization of the eyeglasses, necklaces, earrings, bracelets, electronic products, watches, false teeth, printing head, razors, pen head, toothbrushes, coins, badge, nipples, cups, tableware, fruit, etc.



Biolab Scientific Ltd.

Trillium Executive Center, East Tower, 675 Cochrane Dr, Markham, Ontario L3R 0B8, Canada
Email: info@biolabscientific.com | Website: www.biolabscientific.com